COMPEX® 

INSTRUCTION FOR USE

MODE D'EMPLOI - 66

ANWEISUNGEN - 135

ISTRUZIONI - 205

INSTRUCCIONES - 272

COMPEX® 

INSTRUCTION FOR USE

TABLE OF CONTENTS

1. Explanation of symbols	02
2. How does electrostimulation work?	03
3. Device function	05
Contents of kits and accessories	05
Device description	06
Battery insertion	07
Connections	07
Preliminary settings	08
Choosing a category	09
Selecting a programme	09
Adjusting stimulation intensities	11
Programme progression	12
End of a programme	13
Battery level and charging	13
4. Troubleshooting	14
5. Device maintenance	16
6. Technical specifications	17
7. EMC Table	19
8. Programmes and specific applications	23
9. Table with placements of electrodes and stimulation positions	61

It is strongly recommended to read these instructions and the contra-indications and safety measures carefully before using your stimulator.

1. EXPLANATION OF SYMBOLS

 	See the instructions
	Type BF applied parts.
	Manufacturer's name and address and date of manufacture
	This device must be separated from household waste and sent to special collection facilities for recycling and recovery
	The stand-by button is multi-functional
	Protect from sunlight
	Store in a dry place
IP20 on the unit	This is an indication for protection against ingress of water and particulate matter. The mark IP20 on your unit means: your unit is protected against solid foreign objects of 12.5mm dia and greater. Not protected against water
IP02 on the case	IP02 on the carrying case means: Protected from the ingress of water droplets from a shower of rain.
	Latex-free
REF	Reference number
LOT	Batch number

2. HOW DOES ELECTROSTIMULATION WORK?

Electrostimulation involves stimulating nerve fibres by electrical impulses transmitted by electrodes. The electrical impulses produced by Complex stimulators are high-quality impulses that are safe, comfortable and effective and stimulate various types of nerve fibres:

1. Motor nerves to stimulate a muscular response, referred to as electrical muscle stimulation (EMS).
2. Certain types of sensitive nerve fibres to obtain analgesic effects or pain relief.

1. STIMULATION OF MOTOR NERVES (EMS)

With voluntary activity, the brain orders muscles to contract and a command is then sent to nerve fibres in the form of an electrical signal. This signal is then sent to muscle fibres, which contract. The principle of electrostimulation accurately reproduces the process involved in a voluntary contraction. The stimulator sends an electrical impulse to nerve fibres to excite them. This excitation is then transmitted to muscle fibres and results in a basic mechanical response (= a twitch). This is the basic requirement for muscular contraction. The muscular response is to all intents and purposes identical to the muscular work controlled by the brain. In other words, the muscle does not distinguish between a command sent by the brain or the stimulator.

Programme settings (number of impulses per second, duration of contraction, rest time, total programme duration) subject the muscle to various types of work, depending on the muscle fibre. Various types of muscle fibres can be identified according to their respective contraction speeds: slow, intermediate and fast fibres. A sprinter clearly has more fast fibres and a marathon runner has more slow fibres. With good knowledge of human physiology and full control of stimulation settings of the various programmes, muscular work can be specifically pinpointed to achieve the desired objective (muscular strengthening, increased blood circulation, firming, etc.).

2. STIMULATION OF SENSORY NERVES

Electrical impulses can also excite sensory nerve fibres to obtain analgesic effects or pain relief. Stimulating tactile sensory nerve fibres blocks pain being transmitted to the nervous system. Stimulating another type of sensory fibre increases the production of endorphins and therefore reduces pain. With pain relief programmes, electrostimulation can be used to treat acute or chronic localised pain and muscle pain.

Caution: Do not use pain relief programmes for an extended period without medical advice.

BENEFITS OF ELECTROSTIMULATION

Electrostimulation is a very effective method for making muscles work:

- With a significant improvement in various muscular qualities
- With no cardiovascular or mental fatigue
- With limited stress exerted on joints and tendons. Electrostimulation thereby enables more muscular work than voluntary activity.

For optimum results, Compex recommends supplementing your electrostimulation sessions with other commitments, such as:

- Regular physical exercise
- A balanced and healthy diet
- A balanced lifestyle

3. DEVICE FUNCTION


CONTENTS OF KITS AND ACCESSORIES

	REF	QTY
STIMULATOR	506100	1
CHARGER	6830XX	1
SET OF 4 SNAP CABLES	601131	1
BAG OF SMALL 5X5 ELECTRODES	42215-8	1
BAG OF LARGE 5X10 ELECTRODES	42216-4	1
ELECTRODE PLACEMENTS POSTER	880303	1
QUICK START GUIDE	885216	1
WARNING LEAFLET	M307816	1
CARRY CASE	680029	1
BATTERY PACK	94121X	1

Only use this device with cables, electrodes, battery, power adaptor and accessories recommended by Complex.

DEVICE DESCRIPTION

A Standby button

B Sockets for the 4 stimulation cables

C Stimulation cables

D +/- buttons for the 4 stimulation channels

E I-button, which allows:

- The intensities to be increased on several channels at the same time

F Charger socket (slide the red cover to the right to reveal the charger connector)

G Battery compartment

H Socket for belt clip

BATTERY INSERTION

Open the battery compartment cover and insert the battery, with the label pointing upwards so that the + and - terminals are opposite the device's contacts. Then replace the cover. If the device is not going to be used for over 3 months, ensure that the battery is fully charged. If the device is not going to be used for over 6 months, ensure that the battery is fully charged and remove it from the stimulator. Switch the stimulator off before removing the battery.

CONNECTIONS


CHARGER CONNECTION

Remove all stimulation cables from the stimulator before recharging it. Connect the charger to a wall socket and connect the stimulator by sliding the red cover to the right to reveal the charger connector. It is strongly recommended to fully charge the battery before its first use to improve its performance and life span.

CONNECTING THE CABLES

The cables of the stimulator connect to the 4 sockets at the front of the device.

PRELIMINARY SETTINGS

Before your first use, you can define various settings. This options screen can then be displayed by switching the device off and by pressing and holding the on/off button for at least 2 seconds.

B Press the +/- button of channel 1 to select the language to be used.

C Press the +/- button of channel 2 to adjust screen contrast.

D Use the +/- button of channel 3 to adjust the volume.

E Press the +/- button of channel 4 to adjust the backlighting.

ON: backlighting always active.

OFF: backlighting always inactive.

AUTO: backlighting activated every time a button is pressed.

A Press the on/off button to confirm and save your selections. Settings are immediately applied.

CHOOSING A CATEGORY

N.B. The following screens are generic examples but they work in the same way regardless of the device that you have.

To switch the stimulator on, briefly press the on/off button.

Before selecting a programme, you should select the desired category.

- A** Press the on/off button to switch the device off.
- B** Press the +/- button of channel 1 to select a category.
- E** Press the +/- button of channel 4 to confirm your selection.


SELECTING A PROGRAMME

- A** Press the on/off button to return to the previous screen.
- B** Press the +/- button of channel 1 to select a programme.
- E** Press the +/- button of channel 4 to confirm the selection.

PERSONALISING A PROGRAMME

N.B. The programme personalisation screen is not available for all programmes


- A** Press the on/off button to return to the previous screen.
- B** Press the +/- button of channel 1 to select the muscular group that you want to stimulate.
- C** Press the +/- button of channel 2 to skip the warming-up.
- D** Press the +/- button of channel 3 to select the stimulation cycle.
- E** Press the +/- button of channel 4 to confirm the selection and start the stimulation session.

N.B. It is recommended to start with the 1st cycle and progress to the next level when the cycle is complete, normally after 4 to 6 weeks of stimulation based on 3 sessions per week. It is also important to have reached significant stimulation intensities in sessions before going on to another cycle. At the end of the cycle you can either start a new cycle or do maintenance training based on one session per week.

ADJUSTING STIMULATION INTENSITIES

When you start a programme, you are asked to increase the stimulation intensities. This step is essential for a successful session.

1 Programme duration in minutes and seconds

2 Programme progression bar. For details of how it works, see the following paragraph: "Programme progression"

A Press the on/off button to switch the unit into Pause mode.


B C D E The four channels flash, going from + to 000. The stimulation intensity must be increased to be able to start stimulation. For this, press the + buttons of the relevant channels until the desired setting is reached.

N.B. To increase the intensities on several channels simultaneously, press the I-button and then increase the intensities. The interdependent channels are displayed in white on a black background.

PROGRAMME PROGRESSION

Stimulation starts properly once the stimulation intensity has been increased.

The examples below explain the general rules. Depending on the programme, there may be slight differences.

- 1** Time remaining (in minutes and seconds) until the end of a programme
- 2** The duration bar showing the duration of the contraction and duration of active rest is only shown during the work sequence
- 3** Session sequences
- 4** Warm-up
- 5** Work period
- 6** Relaxation


A Press the on/off button to temporarily interrupt the programme. To resume, simply press the +/- button of channel 4. The session will resume at 80% of the intensity level being used before it was interrupted.

B C D E The various intensities reached during the contraction phase are shown by a series of black vertical bars; the intensities of the rest phase are illustrated by hatched bars.

Please note that the stimulation intensities for the active rest phase are automatically set to 50% of the contraction intensities. They can be changed during the rest phase. Once changed, they are completely independent of the contraction intensities.

END OF A PROGRAMME

At the end of a session, the following screen is displayed. To stop the stimulator, press the on/off button.

BATTERY LEVEL AND CHARGING

Battery performance depends on the programme and the stimulation intensity applied. It is strongly recommended to fully charge the battery before its first use to improve its performance and life span. Always use the charger supplied by Complex to recharge the battery.

If the device is not going to be used for over 3 months, ensure that the battery is fully charged. If the device is not going to be used for over 6 months, ensure that the battery is fully charged and remove it from the stimulator. Switch the stimulator off before removing the battery.

BATTERY LEVEL

The battery charge level is shown by a battery icon in the bottom left of the screen. The battery icon flashes when the battery is completely flat. The device can then no longer be used. Recharge it immediately.

RECHARGING

Remove all stimulation cables from the stimulator before recharging it. Connect the charger to a wall socket and connect the stimulator by sliding the red cover to the right to reveal the charger connector.


The charging menu shown below appears automatically.

The charging duration is displayed on the screen. As soon as charging is complete, the battery flashes. Disconnect the charger: the stimulator will switch off automatically.

4. TROUBLESHOOTING

ELECTRODE FAULT

The device emits a tone and alternately displays the pair of electrodes symbol and an arrow pointing towards the channel on which a problem has been detected. In the above example, the stimulator has detected an error on channel 1.

Check that the electrodes are connected to this channel.

If electrodes are old, worn and/or if the contact is poor: try using new electrodes.

Try using the stimulation cable on a different channel. If the cable is still faulty, replace it: www.Complex.com

STIMULATION DOES NOT PRODUCE THE USUAL SENSATION

Check that all settings are correct and check electrodes are properly positioned.

Change the position of the electrodes slightly.

STIMULATION CAUSES DISCOMFORT

Electrodes lose their adhesive capacity and no longer provide suitable contact with skin.

Electrodes are worn and must be replaced.

Change the position of the electrodes slightly.

5. DEVICE MAINTENANCE

GUARANTEE

See enclosed leaflet.

MAINTENANCE

Your stimulator does not require neither calibration nor periodic maintenance. However when needed clean with a soft cloth and solvent-free alcohol-based detergent. Use as little liquid as possible to clean the device. Do not disassemble the stimulator or the charger as they contain high-voltage components which could cause electrocution. This must be carried out by Complex-approved technicians or repair services. If your stimulator contains parts that appear to be worn or faulty, please contact the closest Complex customer service centre.

CONDITIONS FOR STORAGE/TRANSPORT AND USE

	STORAGE AND TRANSPORT	USE
TEMPERATURE	-20° C to 45° C	0° C to 40° C
MAXIMUM RELATIVE HUMIDITY	75%	30% to 75%
ATMOSPHERIC PRESSURE	from 700 hPa to 1060 hPa	from 700 hPa to 1060 hPa

Do not use in areas at risk of explosion.

DISPOSAL

Batteries must be disposed of in accordance with national regulatory requirements in force. Any product bearing the WEEE label (a crossed-out wheeled bin) must be separated from household waste and sent to special collection facilities for recycling and recovery.

6. TECHNICAL SPECIFICATIONS

GENERAL INFORMATION

94121x rechargeable Nickel-metal hydride (NiMH) battery (4.8 V / \geq 1200 mA/h).

Battery chargers: only battery chargers with the part number 6830xx can be used to recharge the batteries supplied with the stimulator.

Product and accessories (excluded battery) expected service life: 5 years

Electrode shelf life: refer to electrodes bag

NEUROSTIMULATION

All the electrical specifications are supplied with an impedance from 500 to 1000 ohms per channel.

Channels: four independent and individually adjustable channels, electrically isolated from each other.

Impulsion form: constant rectangular current with pulse compensation to eliminate any direct current component to prevent residual polarisation at skin level.

Maximum pulse intensity: 120 mA.

Pulse intensity increments: manual adjustment of stimulation intensity from 0 to 999 (energy) by minimum increments of 0.5 mA.

Pulse amplitude: from 60 to 400 μ s.

Maximum electrical charge per pulse: 96 microcoulombs (2 x 48 μ C, compensated).

Standard pulse ramp-up time: 3 μ s (20 %-80 % of maximum current).

Pulse frequency: 1 to 150 Hz.

INFORMATION ABOUT ELECTROMAGNETIC COMPATIBILITY (EMC)

The stimulator is designed to be used in typical domestic environments approved in accordance with the safety standard EMC EN 60601-1-2.

This device emits very low levels in the radio frequency (RF) interval and is therefore not likely to cause interference with nearby electronic equipment (radios, computers, telephones, etc.).

The stimulator is designed to withstand foreseeable disturbances originating from electrostatic discharge, magnetic fields from the power supply or radio frequency emitters.

However, it is not possible to guarantee that the stimulator will not be affected by powerful RF (radio frequency) fields originating, for example, from mobile phones.

For more detailed information about electromagnetic emissions and immunity, please contact Complex.

STANDARDS

To ensure your safety, the stimulator has been designed, manufactured and distributed in accordance with the requirements of the amended European Directive 93/42/EEC covering medical devices.

The stimulator also complies with the IEC 60601-1 standard covering general safety requirements for electromedical devices, with the IEC 60601-1-2 standard covering electromagnetic compatibility and the IEC 60601-2-10 standard covering special safety requirements for nerve and muscle stimulators.

In accordance with current international standards in force, a warning must be given about applying electrodes to the thorax (increased risk of cardiac fibrillation).

The stimulator also complies with Directive 2012/19/EU covering waste electrical and electronic equipment (WEEE).

7. EMC TABLE

The Complex stimulator needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided.

Portable and mobile RF communications equipment can affect the Complex stimulator.

The use of Accessories, transducers, and cables other than those specified by the manufacturer, may result in increased Emissions or decreased Immunity of the Complex stimulator.

The Complex stimulator should not be used adjacent to or stacked with other equipment and if adjacent or stacked use is necessary, the Complex stimulator should be observed to verify normal operation in the configuration in which it will be used.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS		
The Complex stimulator is intended for use in the electromagnetic environment specified below. The customer or the user of the Complex stimulator should assure that it is used in such an environment.		
EMISSION TESTS	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
Radiated Emissions CISPR 11	Group 1	The Complex uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted Emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The Complex is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations IEC 61000-3-3	Complies	

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The Complex stimulator is intended for use in the electromagnetic environment specified below. The customer or the user of the Complex stimulator should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV (power lines) Not Applicable (I/O lines)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Line to Line (or Neutral) ±2kV line(s) to earth	±1kV Line to Line Not Applicable (Line to Earth)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 a	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Complex stimulator requires continued operation during power mains interruptions, it is recommended that the Complex stimulator be powered from an uninterrupted power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

a UT is the a.c mains voltage prior to application of the test level.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The Complex stimulator is intended for use in the electromagnetic environment specified below. The customer or the user of the Complex stimulator should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Complex stimulator, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ <p>$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz</p>
Radiated RF IEC 61000-4-3	3 V/m 1.4 GHz to 2,7 GHz 10 V/m 26MHz to 1GHz	3V/m 10V/m	<p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Complex stimulator is used exceeds the applicable RF compliance level above, the Complex stimulator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Complex stimulator.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE COMPEX STIMULATOR

The CompeX stimulator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CompeX stimulator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CompeX stimulator as recommended below, according to the maximum output power of the communications equipment.

RATED MAXIMUM OUTPUT POWER OF TRANSMITTER W	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER M		
	150 KHZ TO 80 MHZ $D = 1.2 \sqrt{P}$	80 MHZ TO 800 MHZ $D = 1.2 \sqrt{P}$	800 MHZ TO 2,5 GHZ $D = 2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

8. PROGRAMMES AND SPECIFIC APPLICATIONS

The applications that follow are given by way of example. They provide a better understanding of how electrostimulation sessions can be combined with voluntary activity. These protocols will help you to determine the best procedure to follow depending on your needs (choice of the programme, muscular group, duration, placement of the electrodes, body position).

You can of course choose a muscular group different from the one indicated, depending on the activities you practice or the body area you wish to stimulate. The specific applications provide you information on the placement of electrodes and the stimulation position to adopt. This information is represented by numbers (placement of electrodes) and letters (stimulation position). They refer to the “Table with placements of electrodes and stimulation positions” in chapter 9 and to the poster with the pictures for the placements of electrodes and the pictograms for the body positions. The table also indicates the stimulation intensity to use, according to the programme.

SPORT CATEGORY

Neglected for many years, muscle preparation has today become indispensable for the competitive athlete. In this respect, muscular electrostimulation is a complementary training technique widely used by an increasing number of athletes aiming to improve their level of performance. Increasing the maximum strength of a muscle, developing muscular volume, increasing the explosive strength of muscles or improving the capacity of muscle fibres to sustain effort over long periods of time are objectives that differ according to the sporting discipline being practised. Ensuring optimal muscle preparation immediately before competition, combining electrostimulation with voluntary muscle training, optimising the effects of training techniques such as stretching, reproducing the muscular stress resulting from “plyometric” training or imposing a “restoration” activity on muscles is easily accessible today thanks to the high specificity of the new programmes offered by your Complex. The programmes of the Complex Sport category are designed to improve the muscle qualities of competitive athletes training seriously on a voluntary basis for at least five hours a week. Use of programmes of the Sport category is not suitable for atrophied muscles that have suffered any kind of pathological process. For such muscles, it is essential to use programmes of the Rehabilitation category (see section “Rehabilitation category” of this manual). The protocols that follow are given by way of example only. They provide a better understanding of how electrostimulation sessions can be combined with voluntary training.

At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached

Programmes	Effects	Uses
Potentiation	<ul style="list-style-type: none"> • Increased speed of contraction and power gain • Less nervous effort to attain maximum strength 	<ul style="list-style-type: none"> • For optimum muscle preparation immediately before competition
Endurance	<ul style="list-style-type: none"> • Improved absorption of oxygen by the stimulated muscles • Improved performance for endurance sports 	<ul style="list-style-type: none"> • For athletes who wish to improve their performance in long-duration sporting events
Resistance	<ul style="list-style-type: none"> • Improved absorption of oxygen by the stimulated muscles • Improved performance for endurance sports 	<ul style="list-style-type: none"> • For athletes who wish to improve their performance in long-duration sporting events
Strength	<ul style="list-style-type: none"> • Increased maximum strength • Increased rate of muscular contraction 	<ul style="list-style-type: none"> • For competitive athletes practising a discipline that requires strength and speed
Explosive strength	<ul style="list-style-type: none"> • Increases the speed at which the level of strength is attained • Improved efficiency of explosive actions (jumping, sprinting, etc.) 	<ul style="list-style-type: none"> • For athletes practising a discipline in which explosive strength is an important factor in performance
Hypertrophy	<ul style="list-style-type: none"> • Increased muscle volume • Increased muscle resistance 	<ul style="list-style-type: none"> • For competitive athletes who wish to increase their muscle mass • For body-builders

Programmes	Effects	Uses
Fartlek	<ul style="list-style-type: none"> • Training and preparation of muscles for all kinds of muscular work (endurance, resistance, strength, explosive strength) through different working sequences 	<ul style="list-style-type: none"> • At the beginning of the season to “re-set” the muscles after a rest period and before more intensive and specific training • During the season for those who do not have any preferences for a single kind of performance and prefer to exercise their muscles with different work-out methods
Stretching	<ul style="list-style-type: none"> • To optimise the effects of the voluntary technique of stretching, by means of a reduction of muscular tonus obtained by specific activation of the antagonist of the stretched muscle (reciprocal inhibition reflex) 	<ul style="list-style-type: none"> • For all athletes who wish to maintain or improve their muscular elasticity • To be used after all training sessions or during a specific stretching session
Active recovery	<ul style="list-style-type: none"> • Marked increase in blood flow • Accelerates the elimination of waste from the muscular contraction • Endorphinic effect (see Pain category) • Relaxing effect 	<ul style="list-style-type: none"> • To improve and accelerate the muscle recovery after an intensive exercise • To be used during the first 3 hours after each session of intensive training or after a competition
Regeneration	<ul style="list-style-type: none"> • Analgesic effect through the release of endorphins • Marked increase in blood flow encouraging oxygenation and drainage • Activation of the oxydative metabolism • Reactivation of proprioceptive pathways 	<ul style="list-style-type: none"> • To be used the day after competition as recovery training or as a complement to this type of training, the intensity of which can then be reduced

SPECIFIC APPLICATIONS**Use of Potentiation programme to optimise the effects of explosive strength (sprints, jumps, squash, football, basketball, etc.) immediately before specific training or competition**

The Potentiation programme should be used on the key muscles involved in the discipline practised. In this example, the key muscles of the sprinter (quadriceps) will be stimulated. For other disciplines the choice of muscles to be stimulated may be different.

The Potentiation programme must not take the place of the voluntary warm-up usually performed before the competition. Activation of the cardiovascular system, short accelerations becoming progressively faster, practice starts and stretching will therefore be carried out by the athlete according to his normal routine. A Potentiation session of short duration (approximately 3 minutes) will be applied on the sprinter's quadriceps immediately before the start of his race (or races, in the event of qualifying events). The specific muscular activation of the Potentiation programme allows the maximum level of performance to be attained in the very first seconds of the race.

Cycle duration: 1x
Always carry out the session as soon as possible before the start and at all events not more than 10 minutes before the start. After 10 minutes, the phenomenon of potentiation of muscle fibres quickly wears off
Programme: Potentiation 8G

Preparation for a cyclist training three times a week who wishes to improve endurance performance

Effort expended over a long period of time makes demands on the aerobic metabolism, for which the decisive factor is the quantity of oxygen consumed by the muscles. To progress in endurance, it is therefore necessary to increase as much as possible the supply of oxygen to the muscles stimulated by this type of effort. Because oxygen is conveyed by the blood, it is essential to have an efficient cardiovascular system, due to voluntary training under certain conditions. However, a muscle capacity to consume the oxygen it receives (oxydative capacity) can also be improved by following a specific work regime.

The Endurance programme of the Sport category leads to a significant improvement in the consumption of oxygen by muscles. Combining this programme with the Capillarization programme (Vascular category), which develops the network of intramuscular capillaries, is particularly beneficial and allows endurance athletes to improve their performance levels.

Cycle duration: 8 weeks, 5 x/week.
<i>Eg for 1 week</i>
Mon: Rest
Tue: 1 x Endurance 8G
Wed: Cycle training 1h30 (moderate speed), then 1x Capillarization 8◇
Thu: 1 x Endurance 8G
Fri: Rest
Sat: Cycle training 60' (moderate speed), then 1 x Endurance 8G
Sun: Cycling 2h30 (moderate speed), then 1 x Capillarization 8◇
Programmes: Endurance 8G and Capillarization 8◇

Preparation for a runner training three times a week who wishes to progress in endurance (half-marathon, marathon)

Running as many miles as possible is essential in order to improve performance in endurance events. However, the strain on tendons and joints this type of training causes is today universally recognised. Integrating Compex electrostimulation into the training of the long-distance runner offers an excellent alternative that can help to overcome this problem. The Endurance programme, which improves the muscles' ability to absorb oxygen, and the Capillarization programme (Vascular category), which develops the capillaries in the muscles, allow greater endurance to be achieved while limiting the weekly mileage and therefore the risk of injury.

Cycle duration: 8 weeks, 5 x/week
<i>Eg for 1 week</i>
Mon: Rest
Tue: 1 x Endurance 8G
Wed:
- Voluntary warm up training 20', then 1-2 series of 6 x [30" fast / 30" slow]
- Slow jogging 10' at the end of the session, then 1 x Capillarization 8◇
Thu: 1 x Endurance 8G
Fri: Rest
Sat: Loose jogging 60', then 1 x Endurance 8G
Sun: Extended run 1h30 (moderate speed), then 1 x Capillarization 8◇
Programmes: Endurance 8G and Capillarization 8◇

Pre-season preparation of lactic capacity for a resistance sport with three active training sessions per week (800 metres, track cycling, etc.)

Example of planning to develop the lactic capacity (resistance) of the quadriceps. For other disciplines, the choice of muscles to be stimulated may be different. During pre-season preparation for sports that make great demands on the lactic anaerobic system (intense effort sustained over as long as possible), it is essential not to neglect specific muscle preparation. Stimulation of the quadriceps (or a different key muscle group depending on the discipline practised) by means of the Resistance programme results in improved anaerobic power, as well as greater muscle tolerance to high concentrations of lactates. The practical benefits are evident: improved performance thanks to better muscular resistance to fatigue for exercises of the lactic anaerobic type. In order to optimise the effects of this preparation, you are advised to supplement it with Capillarization sessions, carried out during the week leading up to the competition (see this part: "Programmes and specific applications", "Vascular category", "Achieving peak form before competition for resistance sports (eg: 800 metres, 1500 metres, etc.)"). An Active recovery session is recommended after the most intensive training. It accelerates the rate of muscle recovery and reduces fatigue during the period of the season when the amount of training is high.

Cycle duration: 6-8 weeks, 4 x/week
<i>Eg for 1 week</i>
Mon: 1 x Resistance 8G
Tue: Voluntary training on the track
Wed: 1 x Resistance 8G
Thu: Voluntary intensive training on the track, then 1 x Active recovery 8◇
Fri: Rest
Sat: Rest
Sun: 1 x Resistance 8G, followed by session voluntary training on the track
Programmes: Resistance 8G and Active recovery 8◇

Preparation for a cyclist training three times a week who wishes to improve his power

Developing the strength of the thigh muscles is always beneficial for the competitive cyclist. Certain forms of training on the bicycle (hill work) can make a contribution in this respect. However, results will be more spectacular if additional training using Compex muscular stimulation is undertaken at the same time. The special regime of muscle contractions of the Strength programme and the large amount of work to which the muscles are subjected allow a significant increase of the strength of the thigh muscles. Moreover, the Active recovery programme, carried out within three hours of the most intensive training sessions, encourages muscular recuperation and makes it possible to follow on with qualitative training under optimum conditions.

Cycle duration: 8 weeks, 5 x/week
<i>Eg for 1 week</i>
Mon: Rest
Tue: 1 x Strength 8G
Wed:
- Cycle training 45' (moderate speed), then 5-10 times on a 500-700 m hill (rapidly)
- Recovery during descent
- Inactivity 15-20', then 1 x Active recovery 8◇
Thu: 1 x Strength 8G
Fri: Rest
Sat: Cycle training 60' (moderate speed), then 1 x Strength 8G
Sun:
- Cycling 2h30-3 h (moderate speed)
- Muscle strengthening on hills (use of a high gear ratio remaining seated), then 1 x Active recovery 8◇
Programmes: Strength 8G and Active recovery 8◇

Preparation for a swimmer training three times a week who wishes to improve his swimming power

In swimming, developing the propulsive force of the upper limbs is an important factor in improving performance. Certain forms of voluntary training practised in the water can contribute to this. However, integrating Compex muscular stimulation into the voluntary training programme makes it possible to achieve far better results. The special muscular contraction regime of the Strength programme and the large amount of work to which muscles are subjected will allow you to increase significantly the strength of the latissimus dorsi, key muscles for the swimmer. Moreover, the Active recovery programme, carried out within three hours of the most intensive training, encourages muscular recuperation and makes it possible to follow on with qualitative training under optimum conditions.

Cycle duration: 8 weeks, 5 x/week
<i>Eg for 1 week</i>
Mon: Rest
Tue: 1 x Strength 18C
Wed:
- 20-30' swimming training (different styles), then 5-10 times 100 m with pull-boy
- Recovery 100 m backstroke
- Inactivity 15', then 1 x Active recovery 18◇
Thu: 1 x Strength 18C
Fri: Rest
Sat: 1 h swimming training including some technical work, then 1 x Strength 18C
Sun:
- 20-30' swimming training (different styles), then 5-10 times 100 m with paddles
- Recovery 100 m backstroke
- Inactivity 15', then 1 x Active recovery 18◇
Programmes: Strength 18C and Active recovery 18◇

Pre-season preparation for a team sport (football, rugby, handball, volleyball, etc.)

Planning example to develop the strength of the quadriceps. Depending on the sport being practised, a different muscle group may be chosen. During the preparatory pre-season period for team sports, it is essential not to neglect specific muscle preparation. In most team sports, the qualities of speed and strength make all the difference. Stimulation of the quadriceps (or a different key muscle depending on the sport concerned) by means of the Compex Strength programme will result in increased speed of contraction and muscular strength. The practical benefits will be obvious: improved speed off the mark and during movement, jumping, shooting power, etc. An Active recovery session, carried out after the most intensive training, accelerates muscular recuperation and reduces fatigue that has built up during the season when the workload is high.

Cycle duration: 6-8 weeks, 6 x/week
<i>Eg for 1 week</i>
Mon: 1 x Strength 8G
Tue: Collective training, then 1 x Active recovery 8◇
Wed: 1 x Strength 8G
Thu: Collective training, then 1 x Active recovery 8◇
Fri: 1 x Strength 8G
Sat: Rest
Sun: Collective training or friendly match, then 1 x Active recovery 8◇
Programmes: Strength 8G and Active recovery 8◇

Maintaining results achieved during preparation for team sports during the competitive period (football, rugby, handball, volleyball, etc.)

This example concerns only athletes who have completed a full cycle of training by electrostimulation (at least 6 weeks) as part of their pre-season preparation. The weekly session of stimulation with the Strength programme should be carried out on the same muscle groups as those stimulated during the preparation period (in our example, the quadriceps). During the season, when matches are played regularly, care should be taken not to overtrain the specific musculature. On the other hand, nor should the benefits of preparation be lost by suspending stimulation training for too long a period. During this period of competition, muscle qualities should be maintained by means of one weekly session of stimulation carried out using the Strength programme. It is also essential to leave a sufficiently long interval between this single weekly stimulation session and the day of the competition (3 days minimum). The Active recovery programme, which must be used during the first three hours after the match, as well as after each session of intensive training, restores muscle equilibrium more quickly.

Cycle duration: During the sport season
<i>Eg for 1 week</i>
Mon: Rest
Tue: Collective training, then 1 x Active recovery 8◇ (if training is intensive)
Wed: 1 x Strength 8G
Thu: Collective training, then 1 x Active recovery 8◇ (if training is intensive)
Fri: Rest
Sat: Rest
Sun: Match, then 1 x Active recovery 8◇ (within the 3 hours that follow the competition)
Programmes: Strength 8G and Active recovery 8◇

Pre-season preparation of the explosive strength of the quadriceps for an athlete training three times a week (long jump or high jump, sprinting, etc.)

For other disciplines the choice of muscles to be stimulated may be different. For all sports where the essential performance factor is explosive muscle strength, specific muscle preparation is the main factor in pre-season preparation. The explosive strength of muscles may be defined as the capacity of a muscle to attain a high level of maximum strength as quickly as possible. To develop this quality, voluntary training relies on tiring muscle training sessions that often include the risk of injury, since they are necessarily carried out using heavy weights. Integrating the use of the Explosive strength programme lightens the muscle training sessions while at the same time offering greater benefits and more time for technical work.

Cycle duration: 6-8 weeks, 4 x/ week
<i>Eg for 1 week</i>
Mon: 1 x Explosive strength 8G
Tue: Voluntary training in stadium
Wed: 1 x Explosive strength 8G
Thu: Voluntary training including technical work jumping pit
Fri: 1 x Explosive strength 8G
Sat: Rest
Sun: Voluntary training in stadium followed by 1 x Explosive strength 8G
Programme: Explosive strength 8G

Muscle volume for a body-builder

Despite repeated efforts during their voluntary training, many body-builders encounter difficulties in developing certain muscle groups. The specific stimulation imposed on muscles by the Hypertrophy programme significantly increases the volume of the stimulated muscles. In addition, for a similar session time, the Complex Hypertrophy programme provides a greater volume gain than voluntary training. The additional training imposed by this stimulation programme on muscles not sufficiently receptive to traditional training provides a solution for the harmonious development of all muscle groups without recalcitrant areas. To obtain optimum progress, you are advised to: 1) precede the Hypertrophy sessions with short voluntary training focused on strength; for example 3 series of 5 repetitions at 90% of maximum force; 2) carry out a Capillarization session (Vascular category) directly after the Hypertrophy session.

Muscle volume for a body-builder training three times a week

In this example, we assume that the bodybuilder wants to work on his/her biceps, which he/she considers to be his/her weak point. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. biceps then calves).

Cycle duration: 8 weeks, 3 x/ week
<i>Progression in the levels</i>
Week 1: Hypertrophy level 1
Weeks 2-3: Hypertrophy level 2
Weeks 4-5: Hypertrophy level 3
Weeks 6-8: Hypertrophy level 4
<i>Eg for 1 week</i>
Mon: Rest
Tue: Voluntary training focused on the muscles of the lower limbs, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x Hypertrophy 20D followed by 1 x Capillarization 20◇
Wed: Rest
Thu: Voluntary training focused on the muscles of the trunk, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x Hypertrophy 20D followed by 1 x Capillarization 20◇
Fri: Rest
Sat: Voluntary training focused on the muscles of the upper limbs, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x Hypertrophy 20D followed by 1 x Capillarization 20◇
Sun: Rest
Programmes: Hypertrophy 20D and Capillarization 20◇

Muscle volume for a body-builder training six times a week

In this example, we assume that the bodybuilder wants to work on his/her calves that he/she considers to be his/her weak point. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. calves then biceps).

Cycle duration: 12 weeks, 5 x/ week
<i>Progression in the levels</i>
Week 1: Hypertrophy level 1
Weeks 2-3: Hypertrophy level 2
Weeks 4-5: Hypertrophy level 3
Weeks 6-8: Hypertrophy level 4
Weeks 9-12: Hypertrophy level 5
<i>Eg for 1 week</i>
Mon: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then Hypertrophy 4A followed by 1 x Capillarization 4◇
Tue: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then Hypertrophy 4A followed by 1 x Capillarization 4◇
Wed: Rest
Thu: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then Hypertrophy 4A followed by 1 x Capillarization 4◇
Fri: Rest
Sat: Voluntary training focused on the muscles of the upper limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then 1 x Hypertrophy 4A followed by 1 x Capillarization 4◇
Sun: Rest
Programmes: Hypertrophy 4A and Capillarization 4◇

Resuming activity after a break in training: variable muscle preparation

In this example, we assume that the sportsman wants to work above all on the latissimus dorsi. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. latissimus dorsi then quadriceps). The fartlek comes from Scandinavia. It involves diversified training carried on in a natural setting. During training, different work-out routines are alternated and different muscles are stimulated. For example, after a few minutes of slow jogging, the trainee does a series of accelerations followed by some jumps, before resuming a slower jog, and so on. The objective is to work on different muscular qualities without however emphasising any one in particular. This activity is frequently done either at the beginning of the season, for general muscular reactivation, or done regularly by leisure-time sports enthusiasts who do not want to emphasise any particular kind of muscular performance, but rather want to stay fit and reach the right level in all kinds of muscular work. At the beginning of the season or after a significant break in training, the resumption of physical and/or sports activity should be done progressively and become increasingly specific. It is thus normal to perform a few initial sessions with the aim of making the muscles do all kinds of work to prepare them for later training that will be more intensive and more oriented towards a specific kind of performance. Through these eight sequences that automatically follow each other, the Fartlek programme imposes different kinds of work to the stimulated muscles and thus makes them used to all kinds of effort.

Cycle duration: 1-2 week, 4-6 x/ week
Programme: Fartlek 1&C

Preparation for a footballer wishing to optimise the effects of active stretching of the hamstrings muscles

The Stretching programme consists in stimulating the antagonist (i.e. opposite) muscle of the muscle subjected to stretching in order to exploit a well-known physiological mechanism: reciprocal inhibition reflex. This reflex, based on the proprioceptive sensitivity of muscles, consists of a very marked muscular relaxation. This allows more efficient stretching, since it is carried out on a more relaxed muscle. Stimulation is therefore performed on the muscle opposing the muscle being stretched. This stimulation consists of a progressive appearance and disappearance (slowly and lasting a long time) of contractions, with complete rest between contractions. It is during the contraction (which increases with each new level) that the athlete stretches the chosen muscle group by means of a traditional voluntary stretching technique. In this example, the stimulation is carried out on the quadriceps in order to facilitate stretching of the hamstrings muscles during the contraction phases.

Cycle duration: Throughout the season, according to the frequency of voluntary stretching sessions
According to the normal duration of stretching (this depends on schools of physical training and each person's feelings), choose the level that seems most appropriate. The proposed duration of stretching is: 10 seconds for level 1 12 seconds for level 2 14 seconds for level 3 16 seconds for level 4 18 seconds for level 5
Programme: Stretching 8, Adopt the starting position for the active stretching exercise

Use of the Regeneration programme to eliminate muscle fatigue more quickly (cross-country running, football, basketball, tennis, etc.) and restore good muscular sensations more quickly

The Regeneration programme, known also as the “day-after programme” must be used on key muscles for the discipline being practised. In this example, key muscles for the cross-country runner (quadriceps) will be stimulated. For other disciplines, the choice of muscles to be stimulated may be different. Note that this type of programme is particularly beneficial for all sports where competitions are repeated at frequent intervals, during tournaments and cups in various sports. The stimulation session using the Regeneration programme must be carried out the day after a competition to replace or supplement so-called “restoration” training, which can therefore be less intensive. Contrary to the Active recovery programme, which provokes no tetanic contractions and which must be used during the three hours after the competition or intensive training, the Regeneration programme is a form of light training which, in addition to an analgesic effect and an increased blood flow, aims to impose a small degree of anaerobic training and likewise to provoke slight tetanic contractions that are not tiring, making it possible to reactivate the proprioceptive pathways. Energy pathways are also gently stimulated, allowing their metabolic equilibrium to re-establish itself.

The session is composed of 6 stimulation sequences that follow on automatically:

1st sequence: analgesic effect

2nd sequence: marked increase in blood flow

3rd sequence: tetanic contractions to restore muscular sensation

4th sequence: activation of the oxydative metabolism

5th sequence: marked increase in blood flow

6th sequence: relaxing effect

Cycle duration: To be used throughout the season, according to the frequency of competitions

Programme: Regeneration 8G

PAIN RELIEF CATEGORY

Physical pain is an abnormal and unpleasant sensation caused by an injury, a disorder or incorrect functioning of a part of our organism. It is always a signal sent to us by our bodies, one that should not be ignored, and that in all cases requires us to consult a doctor if it does not disappear quickly. The approach to pain adopted by the medical profession has changed considerably in recent years. Treatment of the cause is always fundamental, however the pain as such must be otherwise removed or at least considerably reduced and made bearable for the patient. The means to combat pain have developed greatly, and there is no longer any hesitation today in using powerful analgesics to improve the quality of life of patients. It is within this context that recourse to electrotherapy has developed. Excitation of the nerve fibres of sensation by means of electrical micro-impulses is often the preferred method of combatting pain. Today, such analgesic electrotherapy is widely used, particularly in rehabilitation medicine and in specialist pain treatment centres. The precision of electrical currents allows the analgesic action to be accurately targeted according to the type of pain. It is up to the user to choose the most appropriate programme for the type of pain and to follow practical recommendations in order to obtain maximum effect.

If the pain is great and/or persistent, you are advised to consult a doctor. Only a doctor can make an accurate diagnosis and put in place therapeutic actions designed to encourage the disappearance of the disorder.

Table of the pathologies

Pathologies	Programmes	References
Neuralgia of the upper limb (brachial neuralgia)	Modulated TENS	See application page 42
Chronic muscular pain (polymyalgia)	Endorphinic	See application page 42
Contracture (eg: localized contracture in external side of the calf)	Decontracture	See application page 43
Chronic muscular pain in the back of the neck (cervical pain)	Cervical pain	See application page 43
Muscular pain in the thoracic region (thoracic back pain)	Thoracic back pain	See application page 44
Muscular pain in the low back region (low back pain)	Low back pain	See application page 44
Sharp and recent muscular pains affecting a muscle in the low back (lumbago)	Lumbago	See application page 45
Chronic elbow pain (epicondylitis = tennis elbow)	Epicondylitis	See application page 45

Table of Pain programmes

Programmes	Effects	Uses
Modulated TENS	Blocks transmission of pain by the nervous system	• All acute or chronic localized pain
Endorphinic	<ul style="list-style-type: none"> • Analgesic action through the release of endorphins • Increased blood flow 	• To combat chronic muscular pain
Decontracture	<ul style="list-style-type: none"> • Reduced muscular tension • Relaxing effect 	• To combat recent and localized muscular pain
Cervical pain	<ul style="list-style-type: none"> • Analgesic action through the release of endorphins • Increased arterial flow 	• Analgesic current specifically adapted to pains in the back of the neck
Thoracic back pain	<ul style="list-style-type: none"> • Analgesic action through the release of endorphins • Increased arterial flow 	• Analgesic current specifically adapted to low back pain (lumbar region)
Low back pain	<ul style="list-style-type: none"> • Analgesic action through the release of endorphins • Increased blood flow 	• Analgesic current specifically adapted to persistent low back pain (lumbar region)
Lumbago	<ul style="list-style-type: none"> • Reduced muscular tension • Relaxing effect 	• Analgesic current specifically adapted to sharp and sudden low back pain (lumbar region)
Epicondylitis	• Blocks transmission of pain by the nervous system	• Analgesic current specifically adapted to persistent pain in the elbow

The use of Modulated TENS, Endorphinic and Decontracture programmes should not be prolonged without medical advice.

SPECIFIC APPLICATIONS

Neuralgia of the upper limb (brachial neuralgia)

Some people suffer from arthritis in the joints of the vertebrae at the back of the neck, or from arthritis or peri-arthritis of the shoulder. These situations often give rise to pain that runs down one arm and is known as “brachial neuralgia”. These arm pains that start in the shoulder or the back of the neck can be reduced with the Compex Modulated TENS programme by following the practical recommendations outlined below.

Cycle duration: 1 week, 1 x/day minimum, then adapt according to how the pain develops
According to requirements, the Modulated TENS programme can be repeated a number of times during the same day
Programme: Modulated TENS 35◇

Chronic muscular pain (polymyalgia)

Some people suffer from muscular pain that affects several muscles or parts of different muscles at the same time. The localization of this chronic pain can vary over time. These continual and diffuse muscle pains are the result of chronic contractures in which acids and toxins accumulate, irritating the nerves and causing pain. The Endorphinic programme is particularly effective against such pain since, in addition to its pain-reducing effect, it increases the blood flow to the contracted muscular groups and removes from them accumulations of acids and toxins. Example: localized pain in the biceps. However, the pain may affect other muscular groups. The practical application illustrated below remains valid, but it is then necessary to place the electrodes on the muscular group concerned.

Cycle duration: 4 weeks, 2 x/day, with a 10 minute break between the 2 sessions
You are advised to consult your doctor if no improvement is observed after the first week of use
Programme: Endorphinic 20◇

Contracture (eg: localized contracture in external side of the calf)

After tiring muscle work, intense training or a sporting competition, certain muscles or certain parts of muscles often remain tense and slightly painful. These are known as muscular contractures which should disappear after a few days with rest, good rehydration, a balanced food intake with mineral salts and application of the Decontracture programme. The phenomenon of contracture frequently affects the calf muscles, but can also occur with other muscles. In this case, simply follow the same practical recommendations as below, but place the electrodes on the muscular group concerned.

Cycle duration: 1 week, 1 x/day
You are advised to consult your doctor if no improvement is observed after the first week of use
Programme: Decontracture 24◇

Chronic muscular pain in the back of the neck (cervical pain)

Positions in which the muscles of the back of the neck remain tense for long periods of time, for example working in front of a computer screen, may be responsible for the onset of pain in the back of the neck or on either side of the base of the neck, in the upper back. These pains are due to a contracture of the muscles, of which the prolonged tension crushes the blood vessels and prevents the blood from supplying and oxygenating the muscle fibres. If this phenomenon is prolonged, there is an accumulation of acid and the blood vessels atrophy. The pain then becomes continuous or appears after only a few minutes spent working in an unfavourable position. These chronic pains in the back of the neck can be treated effectively with the Cervical pain programme, which reactivates the circulation, drains accumulations of acid, oxygenates the muscles, develops the capillaries and relaxes the contracted muscles.

Cycle duration: 4 weeks, 2 x/day, with a 10 minute break between the 2 sessions
You are advised to consult your doctor if no improvement is observed after the first week of use
Programme: Cervical pain 15◇

Muscular pain in the thoracic region (thoracic back pain)

Vertebral arthritis and positions in which the muscles of the spinal column remain under tension for long periods of time are often responsible for the onset of pain in the middle of the back, which is accentuated with fatigue. Pressure of the fingers on the muscles on either side of the spinal column can often trigger a sharp pain. These pains are due to a contracture of the muscles, of which the prolonged tension crushes the blood vessels and prevents the blood from supplying and oxygenating the muscle fibres. If the phenomenon is prolonged, there is an accumulation of acid and the blood vessels atrophy. The pain then becomes continuous or can appear after only a few minutes spent working in an unfavourable position. These chronic back pains can be effectively treated with the Thoracic back pain programme which reactivates the circulation, drains accumulations of acid, oxygenates the muscles, develops the capillaries and relaxes the contracted muscles.

Cycle duration: 4 weeks, 2 x/day with a 10 minute break between the 2 sessions
You are advised to consult your doctor if no improvement is observed after the first week of use
Programme: Thoracic back pain 13◇

Muscular pain in the low back region (low back pain)

Low-back pain is the most frequently encountered pain. In a standing position, the entire weight of the trunk is concentrated on the joints between the last vertebrae and the sacrum. The low back region is therefore under particular strain. The discs between the vertebrae are crushed and the low back muscles contracted and painful. There are very many treatments for the relief of low back pain sufferers; among these, the specific currents of the Compex Low back pain programme provide an appreciable improvement and can even resolve the problem if it is essentially muscular in origin.

Cycle duration: 4 week, 2 x/day, with a 10 minute break between the 2 sessions
You are advised to consult your doctor if no improvement is observed after the first week of use
Programme: Low back pain 12◇

Sharp and recent muscular pains affecting a muscle in the low back (lumbago)

During back movement, for example when lifting something, when turning or when standing up straight after bending down, a sudden pain can be triggered in the low back. Those who suffer from this problem present a contracture of the low back muscles and feel a sharp pain in this region; because they cannot stand completely straight, they remain bent over on one side. All of these symptoms indicate what is known as lumbago, which is principally the result of a sharp and intense contracture of the low back muscles (lumbar region). In a situation such as this, it is always necessary to consult a doctor to receive appropriate treatment. In addition to such treatment, the specific Compex Lumbago programme can help effectively to relax the muscles and remove the pain.

Cycle duration: 4 weeks, 3 x/week
You are advised to consult your doctor if no improvement is observed after the first week of use
Programme: Lumbago 33◇

Chronic elbow pain (epicondylitis = tennis elbow)

All tendons of the muscles that allow us to stretch our hand, wrist and fingers end at the small external bone mass of the elbow (epicondyle). Hand and finger movements therefore transmit tensions that are concentrated in the tendon endings at this bone mass. When hand movements are repetitive, as is the case for painters, tennis players or even those who constantly use the mouse of a computer system, small injuries, accompanied by inflammation and pain, develop in the region of the epicondyle. This is known as “epicondylitis”, which is characterised by pain around the external bone mass of the elbow when pressure is applied or when the forearm muscles are contracted. The Compex Epicondylitis programme supplies specific current to combat this type of pain. It acts effectively as a supplement to rest. However, it is necessary to consult your doctor if the pain gets worse or does not disappear quickly after a few sessions.

Cycle duration: 1 week, 2 x/day minimum, then adapt according to how the pain develops
According to requirements, the Epicondylitis programme can be repeated a number of times during the same day.
Programme: Epicondylitis 36◇

VASCULAR CATEGORY

The low frequency current used of the Vascular category of Compex programmes significantly improves blood circulation in the stimulated area. Indeed many people, more particularly women who remain standing for long periods, suffer from circulatory problems. These mainly affect the legs and are caused by stagnation of the blood and the lymph and are manifested by a feeling of “heavy legs”, swelling, or the dilatation of surface veins. The consequences are multiple: fatigue, tension, pain, lack of oxygenation of tissue and the appearance of varicose veins and oedemas. According to the programme used, the muscular twitches are more or less rapid, separate, and adopt different rhythms. The result is a specific action for each programme, and for this reason you are advised to follow closely the indications of the different treatments so that optimum results can be obtained.

If the symptoms are serious and/or persistent, you are advised to consult a doctor. Only a doctor can establish an accurate diagnosis and put in place all therapeutic measures required to cure the disorder.

Programmes	Effects	Uses
Capillarization	<ul style="list-style-type: none"> • Very marked increase in blood flow • Development of capillaries 	<ul style="list-style-type: none"> • In the period before competition, for those participating in endurance or resistance sports • To improve endurance in those whose fitness level is not high
Heavy legs	<ul style="list-style-type: none"> • Increased venous return • Increased supply of oxygen to muscles • Drop in muscular tension • Elimination of tendency to have cramp 	<ul style="list-style-type: none"> • To remove the feeling of heaviness in the legs arising in unusual situations (standing for long periods, heat, hormonal imbalance linked to the menstrual cycle, etc.)
Cramp prevention	<ul style="list-style-type: none"> • Improved circulation to prevent the onset of nocturnal cramp or cramp after strains 	<ul style="list-style-type: none"> • Cramp mainly affects muscles of the lower limbs, in particular the calf muscles
Lymphatic drainage	<ul style="list-style-type: none"> • Deep massage of the stimulated area • Activation of return lymphatic circulation 	<ul style="list-style-type: none"> • To contrast occasional swelling of the feet and ankles • As a supplement to pressure therapy or manual lymphatic drainage <p><i>Do not use without medical advice in case of oedemas that are the side effect of a pathological problem</i></p>

SPECIFIC APPLICATIONS

Preparation for seasonal endurance activity (eg: hiking, bicycle touring)

Sporting physical activities that last a long time (walking, cycling, cross-country skiing, etc.) require strong muscles with a good capillary circulation so that the muscle fibres can be properly oxygenated. When a long duration physical activity is not practised, or not practised regularly enough, the muscles lose their capacity to effectively consume oxygen and the capillary system becomes rarefied. This lack of muscular quality makes exercising uncomfortable and limits the speed of recuperation, and is responsible for many unpleasant side effects such as drowsiness, contractures, cramp and swelling. To restore to muscles their endurance and ensure full activation of the capillary system, the Compex offers a very efficient mode of stimulation.

Cycle duration: 6-8 weeks, before starting the hike, 3x/week, alternating muscular groups

Programmes: Capillarization 8◇ and 25◇
--

Achieving peak form before competition for resistance sports (eg: 800 metres, 1500 metres, mountain biking, mountain stage in cycling, cycling pursuit on track, 200 metres swimming, slalom)

Sports that require maximum effort lasting between 30 seconds and 5 minutes are the so-called resistance sports. Rapid fibres must work at a capacity close to their maximum and be capable of maintaining this intensive work throughout the event without weakening, in other words the rapid fibres must be resistant. The Capillarization programme, which produces a very significant increase in the blood flow in the muscles, brings about development of the intramuscular capillary system (capillarization). This growth of the capillaries works best if it occurs around the rapid fibres. In this way, the latter's exchange surface with the blood increases, allowing an improved supply of glucose, better diffusion of oxygen and quicker evacuation of lactic acid. Capillarization therefore allows the rapid fibres to be more resistant and to maintain their optimum capacity over a longer period. However, prolonged or too frequent use of this programme may induce a modification of rapid fibres into slow fibres, thereby risking a decline in performance for strength and speed sports. It is therefore important to follow closely the recommendations below concerning the specific application in order to benefit from the positive effects of this treatment.

Cycle duration: 1 week before competition, 2 x/day, with a 10 minute break between the 2 sessions

Programme: Capillarization 8◇

Prevention of contractures in the muscles of the back of the neck for cyclists

During voluntary training or competitions, some athletes experience the problem of contractures in support muscles (for example, the back of the neck for cyclists). Use of the Capillarization programme reduces, and may resolve, this problem. The marked increase in blood circulation and development of the capillary system improve the oxygenation of fibres and their exchanges with the blood. The muscle therefore becomes less prone to contractures.

Cycle duration: 3 weeks, 1 x/day

Programme: Capillarization 15◇

Prevention of the feeling of heaviness in the legs

The feeling of heaviness in the legs is due to a temporary deficiency in the venous return, and not to major organic lesions. Stagnation of the blood mass in the legs is encouraged in certain situations: prolonged periods in a standing position, long periods continuously seated, intense heat, etc. The insufficient supply of oxygen to the tissue, and in particular the muscles, that results from this situation, causes the feeling of heaviness and discomfort in the legs. Use of the Heavy legs programme accelerates the return of blood to the veins and has an important relaxing effect on the painful muscles.

Cycle duration: Use this treatment whenever an occasional feeling of heaviness in the legs occurs

Programme: Heavy legs 25 JL

Prevention of cramp in the calf muscles

Many people suffer from cramp in the calf muscles, which can appear spontaneously during rest at night or as a result of prolonged muscular effort. The phenomenon of cramp may in part be caused by a disequilibrium of blood flow in the muscles. To improve the blood circulation and prevent the onset of cramp, the Compex offers a specific stimulation programme. Use of this programme in accordance with the procedure below will allow you to obtain positive results and limit the occurrence of cramp.

Cycle duration: 5 weeks, 1 x/day, at the end of the day or at night

Programme: Cramp prevention 25◇

Prevention of occasional swelling of the feet and ankles

Insufficient blood circulation in the veins in unusual circumstances frequently gives rise to an accumulation of blood and lymph in the extremities of the lower limbs. This phenomenon, which results in swelling of the ankles and feet, brings with it an unpleasant sense of heaviness and tension in the areas concerned. The deep massage caused by the Lymphatic drainage programme will effectively activate lymph return circulation and thus encourage the elimination of this disorder.

Cycle duration: 8 weeks, 3 x/week
Use this treatment whenever you experience occasional swelling of the feet and ankles
Programme: Lymphatic drainage 25◇

MASSAGE CATEGORY

The programmes of the Compex Massage category subject the muscles of the stimulated region to moderate activity, which produces beneficial effects that help to improve physical comfort and well-being. Remaining in the same working position for a long time (for example, sitting in front of a computer screen), conditions of stress, repeated jostling, insufficient muscular conditioning before physical activity are all very frequent situations that are often responsible for uncomfortable bodily sensations. The automatic progression of the various stimulation sequences, specific to each programme of the Massage category, allows each user to adapt the effects of stimulation to his own needs and objectives.

If strong pain persists, the programmes of the Massage category should never be used for long periods without first seeking medical advice.

Programmes	Effects	Uses
Relaxing massage	<ul style="list-style-type: none"> • Decreased muscular tension • Drainage of the toxins responsible for the exaggerated increase of muscular tone • Effect of well-being and relaxation 	To eliminate uncomfortable or painful sensations, following an exaggerated increase of muscular tone
Reviving massage	<ul style="list-style-type: none"> • Considerable increase of blood flow in the stimulated region • Improvement of tissue oxygenation • Elimination of free radicals 	To contrast effectively sensations of fatigue and localised heaviness
Toning massage	<ul style="list-style-type: none"> • Activation of blood circulation • Recovery of muscular contractile properties • Invigorating effect 	To prepare the muscles in an ideal manner before an unusual/one-time physical activity

SPECIFIC APPLICATIONS

Treatment of uncomfortable muscular tensions in the back of the neck

Remaining in a seated position for long periods, associated with repetitive movements of the upper limbs (as is frequently the case in front of a computer screen), is often responsible for an uncomfortable – or even painful – increase of muscle tension of the back of the neck muscles. Any other condition of stress may also create situations leading to a state of excessive muscular tension, which is often responsible for painful or uncomfortable sensations. The in-depth effects produced by the Relaxing massage programme enable to contrast effectively such painful sensations, with a particularly significant result in terms of relaxation.

Cycle duration: To be used on the back of the neck muscles, whenever an occasional sensation of painful muscular tension is felt; to be repeated, if necessary, in cases of particularly strong muscular tension
--

Programme: Relaxing massage 15◇

Treatment of a localised sensation of heaviness or an occasional state of fatigue

The many stresses of everyday life are often responsible for uncomfortable, or even painful, physical sensations. A circulation slowdown is often the result of insufficient physical activity, frequently aggravated by the need in one's working life to remain in the same position for many hours (in a seated position, for example). Although not serious, this simple "vascular slowdown" is, however, frequently the cause of unpleasant sensations (for example: a sensation of heaviness, often localised in the lower limbs, but also, at times, in another region of the body). The Reviving massage programme produces, in the greatest comfort, a reactivation of the blood circulation that allows an acceleration of tissue oxygenation and the elimination of painful sensations, due to insufficient physical effort.

Cycle duration: To be used on the calf muscles, whenever an occasional sensation of heaviness is felt; to be repeated, if necessary, in cases of persistent discomfort
--

Programme: Reviving massage 25◇

Muscular and circulation conditioning before physical activity

People who regularly practise sports are well aware of the transition, which is often painful, between a rest activity and an occasional physical effort. The objective of the usual warming-up techniques is to respond to this need to activate progressively the physiological functions involved in one-time physical activity. This physiological need is also very desirable for moderate but unusual physical efforts, as is frequently the case for most of us (treks, biking, jogging, etc.). The Toning massage programme offers an ideal benefit in the form of optimal muscular and circulation conditioning before any type of physical effort. It allows us to avoid the uncomfortable sensations – usually experienced during the first few minutes of unusual physical efforts – and to limit the secondary consequences of insufficient preparation (aches, etc.).

Cycle duration: To be used on the muscles most subject to stress in the physical activity involved (in this example: the quadriceps), in the last thirty minutes preceding physical activity.

Programme: Toning massage 8◇

REHABILITATION CATEGORY

The muscular consequences of any pathological problem must be treated by means of specific programmes. Indeed, improving the qualities of a healthy muscle or restoring the potential of a “convalescent” muscle require appropriate work regimes. In order to “redevelop” a convalescent muscle, it is necessary to use the programmes of the Rehabilitation category. Muscle volume is quickly reduced after trauma to a bone or a joint, particularly if the injury is treated by immobilisation and/or surgical intervention. This muscular atrophy can be much more progressive in the case of a degenerative injury (such as arthrosis, for example), since the reduction of muscle activity then often increases slowly and is superimposed on the developing pathology. The phenomenon of atrophy is accompanied by a reduction of muscle strength; however, these disturbances are the consequence of various alterations at the level of the muscle fibres. In a rehabilitation process, it is essential, initially, to treat the atrophy (i.e. regain normal muscle volume) before attempting to increase the strength of the muscle using the Reinforcement programme. The use of programmes of the Rehabilitation category must under no circumstances replace rehabilitation sessions carried out in the presence of the physiotherapist. Although the restoration of initial muscular qualities is a key factor in the process of rehabilitation, other aspects (joint mobility, vigilance, residual pain, etc.) can only be dealt with effectively by a competent health professional. Some pathologies and some post-operative rehabilitations require special precautions during use; you are therefore always advised to seek the opinion of your doctor or physiotherapist before using any of the programmes of the Rehabilitation category.

Osteosynthesis equipment: The presence of osteosynthesis equipment (metallic equipment in contact with the bone: pins, screws, plates, prostheses, etc.) is not a contra-indication for the use of Complex programmes. The electrical current of the Complex are specially designed to have no harmful effect with regard to osteosynthesis equipment.

Programmes	Effects	Uses
Disuse atrophy	<ul style="list-style-type: none"> • Reactivation of the trophic action of muscle fibres damaged during atrophy 	<ul style="list-style-type: none"> • Treatment to combat any reduction of muscle volume: - as a result of trauma necessitating immobilisation - accompanying degenerative lesions to joints
Muscle growth	<ul style="list-style-type: none"> • Increased diameter and capacity of muscle fibres damaged during under-activity or inactivity imposed by some kind of pathological problem • Restoration of muscle volume 	<ul style="list-style-type: none"> • After a period of using the Disuse atrophy treatment, as soon as the muscle shows signs of a slight recovery of volume or tone • Until the virtually complete restoration of muscle volume
Reinforcement	<ul style="list-style-type: none"> • Increase the strength of a muscle previously atrophied • Increase the strength of a muscle affected by a pathological process 	<ul style="list-style-type: none"> • At the end of rehabilitation, when the muscle has regained normal muscle volume • From the onset of rehabilitation for non-atrophied muscles

SPECIFIC APPLICATIONS

Atrophy of the quadriceps as a result of trauma

The quadriceps is a voluminous muscle located in the anterior part of the thigh. It is the main muscle that allows you to stretch the knee; its role is therefore essential for walking, running, climbing stairs, etc. Any trauma affecting a lower limb therefore results in a wasting of this muscle, the reduction of muscle volume being more or less serious depending on the duration of the period of inactivity. This atrophy is normally spectacular when it occurs after trauma to the knee, particularly if the trauma was treated by means of surgical intervention. The programmes of the Rehabilitation category are specifically designed to treat the deterioration of muscle fibres that results from such a process. The progressiveness of the work imposed by the different recommended programmes is decisive in obtaining optimum results.

Cycle duration: 10 weeks

Weeks 1-2: 1 x/ day Disuse atrophy 8G

Weeks 3-8: 1 x/ day Muscle growth 8G

Weeks 9-10: 1 x/ day Reinforcement 8G

Programmes: Disuse atrophy, Muscle growth and Reinforcement 8G

Atrophy of the gluteus as a result of arthrosis of the hip (coxarthrosis)

The pain and stiffening that result from osteoarthritis of the hip normally lead to under-utilisation of the buttock muscles, causing a reduction of the volume and quality of these muscles. The main effect of this atrophy is to produce instability of the pelvis, which causes limping and accentuates the pain by increasing the pressures supported by the joint. Lesions affecting the cartilage unfortunately remain irreversible. However, a remedy for atrophy of the gluteus makes it possible to improve the stability of the joint and therefore make the hip less painful by allowing it to maintain a satisfactory mechanical function.

Cycle duration: 10 weeks, then maintenance

Weeks 1-2: 1 x/ day Disuse atrophy 9ABL

Weeks 3-8: 1 x/ day Muscle growth 9ABL

Weeks 9-10: 1 x/ day Reinforcement 9ABL

Week 11 and following weeks: 1 x/ week Reinforcement 9ABL

Programmes: Disuse atrophy, Muscle growth and Reinforcement 9ABL

Development of the latissimus dorsi to treat and prevent tendinous shoulder pain (rotator cuff syndromes)

The shoulder is a complex joint enabling us to make broad gestures (for example, raising our arms in the air). During some of these movements, the tendons in the shoulder can rub against or be compressed against bony fragments of the joint. When this phenomenon is repeated, or occurs in certain patients who have an unfavourable anatomic constitution, such damage to the tendons causes them to become inflamed and thicker, resulting in a considerable increase in their degree of compression. The pain often then becomes very severe, preventing all movement of the shoulder; it can even set in at night and cause serious sleep deprivation. Appropriate medical treatment can only be put in place by consulting a doctor. However, electrostimulation of the latissimus dorsi by means of specific Complex programmes can reduce the distress to which tendons are subjected by increasing their freedom of movement around the shoulder joint.

Cycle duration: 6 weeks, then maintenance
Weeks 1-2: 1 x/ day Disuse atrophy 18C
Weeks 3-6: 1 x/ day Muscle growth 18C
Week 7 and following weeks: 1 x/ week Muscle growth 18C
Programmes: Disuse atrophy Muscle growth 18C

Development of the abdominal belt to prevent pain in the lumbar region (low back pain)

Low back pain arises most frequently among subjects who do not have sufficient musculature in the abdominal region. These muscles in fact represent a veritable natural “corset”, the role of which is to protect the low back region from excessive stress of all kinds. This is why it is referred to as the abdominal “belt”. After a bout of lumbago, when the pain has stopped (see Section “Pain category”, page 40), a common recommendation is to improve the efficiency of the abdominal and back muscles to prevent any recurrence. The Disuse atrophy programme imposes a large amount of work on the abdominal muscles, without requiring the harmful or even dangerous positions often adopted when voluntary exercises are carried out incorrectly. An abdominal belt that has greater strength and endurance can then satisfactorily fulfil its role of protecting the low back region.

Cycle duration: 4 weeks, then maintenance
Weeks 1-4: 1 x/ day Disuse atrophy 10I
Weeks 5 and following weeks: 1 x/ week Disuse atrophy 10I
Programme: Disuse atrophy 10I

Development of the low back muscles to prevent pain in the lumbar region (low back pain)

Like the muscles of the abdominal region, the muscles of the low back (lumbar muscles) also play a role in protecting the low back region. Subjects whose low back muscles are inadequately effective are particularly prone to low back pain. Once the pain in the low back has disappeared, patients are often advised to strengthen the low back muscles with a view to preventing the recurrence of painful episodes. However, carrying out voluntary back exercises often presents serious difficulties for patients who suffer from low back pain. That is why electrostimulation of the low back muscles using the Disuse atrophy programme is one of the preferred methods to improve the efficiency of these muscles.

Cycle duration: 4 weeks, then maintenance
Weeks 1-4: 1 x/ day Disuse atrophy 14L
Weeks 5 and following weeks: 1 x/ week Disuse atrophy 14L
Programme: Disuse atrophy 14L

Strengthening of the lateral peroneus muscles after ankle sprain

The purpose of the lateral peroneus muscles is to maintain the stability of the ankle joint and prevent it from rotating inwardly. After a sprain, these muscles lose their reflex-contraction capacity together with much of their strength. Regaining competent lateral peroneal muscles after a sprain is a fundamental step, without which recurrence is very probable. To do their job correctly, the lateral peroneals must be strong enough to prevent the foot twisting inwards, but they must also contract reflexively at the precise moment when the heel tilts inwards. To develop both of these aspects, strength and speed of contraction, you should use the Reinforcement programme, which produces efficient lateral peroneal muscles and therefore helps to prevent recurrence.

Cycle duration: 4 weeks, then maintenance
Weeks 1-4: 1 x/ day Reinforcement 2A
Weeks 5 and following weeks: 1 x/ week Reinforcement 2A, if you take part in a dangerous sporting activity
Programme: Reinforcement 2A

FITNESS CATEGORY

Today there are more fitness enthusiasts than ever before, and their number is rising.

Apart from some rare individuals who have real competitive targets, the vast majority have only one aim: to restore their body to peak physical condition or maintain it at that level. Cardio training therefore alternates with more specific exercises to develop or maintain a good quality musculature. With this in mind, the sought-after aims may differ according to who is doing the training: increasing muscle volume to achieve an imposing stature (body-building) or better muscle endurance to improve physical comfort during sustained efforts. Combined with a voluntary physical activity (aerobic exercises in the fitness facility, footing, cycling, swimming, etc.), which becomes more pleasant and therefore more effective, the programmes of the Fitness category enable users to obtain a toned and harmonious figure.

The treatments of the Fitness category are intended to make healthy muscles work, they are not suitable for atrophied muscles that have suffered from any kind of pathological process. For such muscles, it is necessary to use the programmes of the Rehabilitation category.

Programmes	Effects	Uses
Muscle starter	<ul style="list-style-type: none"> • Improvement of the contractile qualities of insufficiently used muscles • Restoration and/or improvement of cellular exchanges at the level of the stimulated muscles 	<ul style="list-style-type: none"> • To reactivate the muscles of sedentary subjects and restore physiological muscular qualities • To complement efficiently a voluntary physical activity practiced to be in good condition
Muscle building	<ul style="list-style-type: none"> • Hypertrophy • Increased diameter of muscle fibres 	<ul style="list-style-type: none"> • For those who wish to gain muscle volume and mass
Aerobic	<ul style="list-style-type: none"> • Increased muscle ability to sustain long-lasting efforts • Improved consumption of oxygen of the stimulated muscles 	<ul style="list-style-type: none"> • To improve physical comfort during aerobic activities • To delay the onset of muscle fatigue during long duration activities • To improve well-being during daily life activities

SPECIFIC APPLICATIONS

Preparing the thighs prior to a week skiing for a sedentary person

Seasonal physical activity, such as skiing, even if practised as a leisure activity, will subject the organism of sedentary subjects to unusual stresses. The thigh muscles (quadriceps) are the muscles that are most used, and these play an essential role in protecting the knee joints. It is for this reason that the pleasure of the first few days' skiing is frequently spoilt by severe muscle sores, and even knee trauma, which can sometimes be serious. The Compex Muscle starter and Muscle building programmes offer ideal muscular preparation for the thigh muscles. You can then go skiing without fear of overdoing things, with enhanced pleasure and safety.

Cycle duration: 6-8 weeks
Week 1-3: 2 x Muscle starter 8G / week
Week 4-6: 3 x Muscle building 8G / week
<i>Progression in the levels</i>
Weeks 1-3: Muscle starter step to the next level every session
Weeks 4-6: Muscle building step to the next level every week
Programmes: Muscle starter and Muscle building 8G

Shaping and developing the shoulders for an active person

Most voluntary physical activities, such as jogging or cycling, do not make intensive use of the shoulder muscles. That is why it is particularly beneficial to offset this under-utilisation by combining Compex sessions with your voluntary training programme. The Muscle building programme imposes a very large quantity of work specifically on muscles of the upper body, bringing about an increase in volume and harmonious development of the shoulders. As opposed to voluntary exercises carried out with heavy weights, that are liable to damage joints and tendons, stimulation with the Compex causes few, if any, strains to joints and tendons

Cycle duration: 4 weeks minimum
<i>Eg for 1 week</i>
Mon: 45'-1h of voluntary physical activity (jogging, swimming, cycling, fitness activity, etc.), then 1 x Muscle building 17H
Tue: Rest
Wed: 1 x Muscle building 18C
Thu: 45'-1h of voluntary physical activity (jogging, swimming, cycling, fitness activity, etc.)
Fri: 1 x Muscle building 17H
Sat: Rest
Sun: 1 x Muscle building 18C
Programmes: Muscle building 18C and 17H

Improve physical comfort and pleasure during long-lasting activity (eg: jogging, biking, leisure swimming)

Many people, even on a non-competitive basis, regularly practice an aerobic type physical activity involving moderate intensity but long-lasting effort. This method is today unanimously considered to be the healthiest approach: it enables the person to improve and maintain cardio-vascular and muscle qualities. Using the Aerobic programme (possibly combined with the Capillarization programme of Vascluar category), these enthusiasts can develop the ability of their muscles to use efficiently the oxygen provided by the organism. More pleasure, a better physical comfort during leisure activities and the supplementary benefit of achieving objectives are thus obtained.

Cycle duration: 8 weeks 4 x/week
<i>Eg for 1 week</i>
Mon: Rest
Tue: 1 x Aerobic 8G
Wed: 45'-1 h of voluntary activity of aerobic type (footing, cycling, swimming, etc.), then according to your possibilities and optionally, 1 x Capillarization 8◇
Thu: 1 x Aerobic 8G
Fri: Rest
Sat: 45'-1 h of voluntary activity of aerobic type (footing, cycling, swimming, etc.)
Sun: If possible, 1 h-1h15 of voluntary activity of aerobic type (footing, cycling, swimming, etc.), then, according to your possibilities and optionally, 1 x Capillarization 8◇
Programmes: Aerobic 8G and Capillarization 8◇

9. TABLE WITH PLACEMENTS OF ELECTRODES AND STIMULATION POSITIONS

The following tables will help you to determine which electrode placement and associated stimulation position to adopt, depending on the chosen programme.

They should be read together with the poster:

- the numbers on the table refer to the electrode placement pictures in the poster;
- the letters refer to the pictograms of the corresponding body positions.

The table also indicates the stimulation intensity to use, according to the programme.

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POSTER


Table with placements of electrodes and stimulation positions

SPORT	Intensity	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
Potentialtion	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	
Endurance	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	
Resistance	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	
Strength	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	
Explosive strength	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	
Hypertrophy	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	
Fartlek	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	
Stretching		Increase progressively the intensities to obtain a tension that is clearly perceived (eg.: 3, 4, 5, 7, 8); Comfortable position																								
Active recovery	*	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇
Regeneration	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	
SPORT	Intensity	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	
Potentialtion	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Endurance	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Resistance	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Strength	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Explosive strength	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Hypertrophy	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Fartlek	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Stretching		Increase progressively the intensities to obtain a tension that is clearly perceived (eg.: 3, 4, 5, 7, 8); Comfortable position																								
Active recovery	*	-	◇	-	-	-	◇	-	-	-	-	-	-	-	-	-	-	-	◇	-	-	-	-	-	-	
Regeneration	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	

PAIN RELIEF	Intensity	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
Modulated TENS	☼	Cover the painful area (eg.: 26, 35, 36, 37, 38, 39); Comfortable position																								
Endorphinic	*	Positive electrode on the painful point (eg.: 12, 13, 15); Comfortable position																								
Decontracture	*	Positive electrode on the painful point (eg.: 24) or standard placement (01-23); Comfortable position																								
Cervical pain	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Thoracic back pain	*	-	-	-	-	-	-	-	-	-	-	-	◇	-	-	-	-	-	-	-	-	-	-	-	-	-
Low back pain	*	-	-	-	-	-	-	-	-	-	-	-	◇	-	-	-	-	-	-	-	-	-	-	-	-	-
Lumbago	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Epicondylitis	☼	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

PAIN RELIEF	Intensity	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	
Modulated TENS	☼	Cover the painful area (eg.: 26, 35, 36, 37, 38, 39); Comfortable position																								
Endorphinic	*	Positive electrode on the painful point (eg.: 12, 13, 15); Comfortable position																								
Decontracture	*	Positive electrode on the painful point (eg.: 24) or standard placement (01-23); Comfortable position																								
Cervical pain	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Thoracic back pain	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Low back pain	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Lumbago	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
Epicondylitis	☼	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

- = Maximum tolerable intensity
- * = Pronounced twitches
- ☼ = Ticking sensation
- ◇ = Comfortable position
- = Non applicable

Letters in the table refer to pictograms for the body positions in the poster

Table with placements of electrodes and stimulation positions

VASCULAR	Intensity	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
Capillarization	*	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇
Heavy legs	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	∟
Cramp prevention	*	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇
Lymphatic drainage	*	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇
VASCULAR	Intensity	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	
Capillarization	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	◇	◇	-	-	-	-	-	
Heavy legs	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
Cramp prevention	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	◇	◇	-	-	-	-	-	
Lymphatic drainage	*	-	-	-	-	-	-	-	-	◇	-	-	-	-	-	-	-	-	◇	◇	-	-	-	-	-	
MASSAGE	Intensity	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
Relaxing massage	*	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	-
Reviving massage	*	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	-
Toning massage	*	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	◇	-
MASSAGE	Intensity	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	
Relaxing massage	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	◇	◇	-	-	-	-	-	
Reviving massage	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	◇	◇	-	-	-	-	-	
Toning massage	*	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	◇	◇	-	-	-	-	-	

Table with placements of electrodes and stimulation positions

REHABILITATION	Intensity	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
Disuse atrophy	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	-
Muscle growth	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	-
Reinforcement	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	-
REHABILITATION	Intensity	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	
Disuse atrophy	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Muscle growth	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Reinforcement	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
FITNESS	Intensity	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
Muscle starter	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	-
Muscle building	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	-
Aerobic	•	◇	A	A	A	IL	E	E	G	ABL	I	I	B	B	B	B	B	H	C	D	D	D	F	F	-	-
FITNESS	Intensity	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	
Muscle starter	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Muscle building	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	
Aerobic	•	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	C	I	-	-	-	-	-	

• = Maximum tolerable intensity
 * = Pronounced twitches
 ☀ = Ticking sensation
 ◇ = Comfortable position
 - = Non applicable
 Letters in the table refer to pictograms for the body positions in the poster

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